



Heater Unit HU 35

System flow chart, version 10


Contents

Heater Unit HU 35.....	4
Country kit.....	4
Water supply tubing.....	6
Fixation.....	6
Hansen couplings.....	6
Cleaning connector.....	6
Disinfectant.....	6
Equipment for descaling/disinfection.....	6
Equipotential bonding cable.....	7
Installation.....	7
Service.....	7
Training.....	7

Heater Unit HU 35, system flowchart, version 10

Heater Unit HU 35

110-115 V




Heater Unit HU 35, 110-115 V
70107.2163

- Including:
 - Coupling kit, straight (70103.3547).
 - Cleaning connector (70102.9306).
- **Country kit must be ordered separately (see below).**
- **Water supply tubing must be ordered separately (see below).**

Note!
Accessories below in green frame are included in Heater Unit HU 35.

230 V




Heater Unit HU 35, 230 V
70107.2162

- Including:
 - Coupling kit, straight (70103.3547).
 - Cleaning connector (70102.9306).
- **Country kit must be ordered separately (see below).**
- **Water supply tubing must be ordered separately (see below).**

Note!
Accessories below in green frame are included in Heater Unit HU 35.

Note! Not all products listed here are available in all countries. For details contact your sales representative.


Note! Please ensure that you select a suitable country kit according to the provided voltage in your country. In case of doubt have an authorized service technician measure the actual supply voltage at the particular power outlet.







Dark grey line = standard system configuration

Country kit

Country Kit, Heater Unit HU 35, 100-115 V




English (EN CAN), Power cable with plug type B		70107.2210	
French (FR CAN), Power cable with plug type B		70107.2213	
Portuguese (BR), Power cable with plug type J		70107.2233	
Spanish (ES – MX, South/ Central America), Power cable with plug type B		70107.2224	









Standard setting of the user interface will be English. It can be changed into German, French, Spanish and Italian. Please see Instructions for Use.









- Including:
 - Mains power inlet cable.
 - Instructions for Use.

Note:
The machine itself is able to handle the whole frequency bandwidth of 50-60Hz within each power supply.

Country kit, Heater Unit HU 35, 230 V


















Chinese (ZH), Power cable with plug type I		70107.2206	
Croatian (HR), Power cable with plug type E+F		70107.2207	
Czech (CS), Power cable with plug type E+F		70107.2208	
Danish (DA), Power cable with plug type K		70107.2209	
Dutch (NL), Power cable with plug type E+F		70107.2226	
English (EN AU), Power cable with plug type I		70107.2227	
English (EN THAI), Power cable with plug type B		70107.2211	
English (EN GB, SING), Power cable with plug type G		70107.2228	

English (EN), Power cable with plug type E+F		70107.2212	
Finnish (FI), Power cable with plug type E+F		70107.2229	
French (FR), Power cable with plug type E+F		70107.2230	
German (DE), Power cable with plug type E+F		70107.2215	
German (DE CH), Power cable with plug type J		70107.2214	
Greek (EL), Power cable with plug type E+F		70107.2216	
Hungarian (HU), Power cable with plug type E+F		70107.2231	
Italian (IT), Power cable with plug type E+F		70107.2217	

110–115 V

Country kit, 230 V

Korean (KO), Power cable with plug type E+F		70107.2218	<input type="checkbox"/>
Norwegian (NO), Power cable with plug type E+F		70107.2334	<input type="checkbox"/>
Polish (PL), Power cable with plug type E+F		70107.2232	<input type="checkbox"/>
Portuguese (PT), Power cable with plug type E+F		70107.2219	<input type="checkbox"/>
Portuguese (BR), Power cable with plug type J		70107.2233	<input type="checkbox"/>
Romanian (RO), Power cable with plug type E+F		70107.2234	<input type="checkbox"/>
Russian (RU, KZ, BY, UA), Power cable with plug type E+F		70107.2220	<input type="checkbox"/>
Slovakian (SK), Power cable with plug type E+F		70107.2235	<input type="checkbox"/>


Slovene (SI), Power cable with plug type E+F		70107.2221	<input type="checkbox"/>
Spanish (ES), Power cable with plug type E+F		70107.2237	<input type="checkbox"/>
Spanish (ES – Central America, PE), Power cable with plug type B		70107.2224	<input type="checkbox"/>
Spanish (ES – South America (AR, UY)), Power cable with plug type I		70107.2222	<input type="checkbox"/>
Spanish (ES – South America (CL)), Power cable with plug type L		70107.2223	<input type="checkbox"/>
Swedish (SV), Power cable with plug type E+F		70107.2225	<input type="checkbox"/>
Turkish (TU), Power cable with plug type E+F		70107.2236	<input type="checkbox"/>

Standard setting of the user interface will be English. It can be changed into German, French, Spanish and Italian. Please see Instructions for Use.

- Including:
 - Mains power inlet cable.
 - Instructions for Use.


Note:
The machine itself is able to handle the whole frequency bandwidth of 50-60 Hz within each power supply.

Water supply tubing



Water supply tubing kit
 1.5 m 70103.3648
 3.0 m 70103.3551
 • 2 water supply tubings, for HU 35 Hansen couplings (self closing, straight or angled).


Fixation



Fixation system for external rail system
 70103.3552

3/6

For SPRINTER CART




Fixation system for Heater Unit HU 35
 70104.8574

For SPRINTER CART XL




Shelf kit for hyperthermia devices
 70105.4144
 • Height adjustable and extensible shelf with railing system.
 • Including four clamps for mounting on the four masts of the SPRINTER CART XL.
 • Suitable for the Heater Unit HU 35.



Tubing holder
 70103.3649


Hansen couplings

This item is included in the Heater Unit HU 35 versions 110-115 V and 230 V.



Hansen coupling kit straight, 1/2"
 70103.3547
 • 2 pieces for in- and outlet (male/female connector).
 To connect the HU 35, via the water supply tubing, to oxygenators.


3/6



Hansen coupling kit angled, 1/2"
 70105.6442
 • 2 pieces for in- and outlet (male/female connector).
 To connect the HU 35, via the water supply tubing, to oxygenators.


Cleaning connector

This item is included in the Heater Unit HU 35 versions 110-115 V and 230 V.



Cleaning connector
 70102.9306
 • Coupling straight, male/male.


Disinfectant



Disinfectant (Clorina)
 70105.6725
 • Powder, 500 g.
 • Active ingredient: Sodium tosylchloramid (Chloramine T).


Equipment for descaling/ disinfection

3/6




Pall-Aquasafe™ disposable water filter
 70107.2030
 • Type AQ31F1R – 31 Day Use.
 • Sterilising grade Supor® membrane 0.2 µm with integrated pre-filtration layer approx 1.0 µm.
 • 2 items per package.


3/6




Quick connect tap adaptors TAP1/2FN
 70107.2031
 • For connecting Pall water filter (1/2 in. female thread) with filling hose.



Filling hose, 3 m length
 70107.2035
 • Outlet Ø 13 mm.
 • With hose barb, 13 mm length, 1/2" external thread and hose clamp.



Equipotential bonding cable

	Equipotential bonding cable	
1.0 m	70105.5328	<input type="checkbox"/>
2.0 m	70105.5348	<input type="checkbox"/>
3.0 m	70105.5349	<input type="checkbox"/>
6.0 m	70105.5350	<input type="checkbox"/>
<ul style="list-style-type: none"> • Connecting leads for potential equalization. • On both ends with assembled right-angled, spring-loaded socket. 		


Installation



Installation

- The equipment should be installed by a Getinge certified technician.
- To be quoted separately.


Service



Service contract


- To be specified by technical dept.

Training



Technical training

- For details contact your local representative.



Application training

- For details contact your local representative.



Getinge is a global provider of innovative solutions for operating rooms, intensive care units, sterilization departments and for life science companies and institutions. Based on our firsthand experience and close partnerships with clinical experts, healthcare professionals and medtech specialists, we are improving the everyday life for people, today and tomorrow.

This document is intended to provide information to an international audience outside of the US

Maquet Cardiopulmonary GmbH · Kehler Str. 31 · 76437 Rastatt · Germany · +49 7222 9320

www.getinge.com